

INSTRUCTION
MANUAL

ANSR: HALO



WELCOME

ANSR: HALO™ is a beautifully simple answer for reducing the signs of aging. ANSR: HALO™ light therapy treatment uses the latest generation of LED technology and specific wavelengths of red and infrared light. These combined wavelengths of light reduce the appearance of fine lines, wrinkles, and discoloration, while helping to even out skin texture and tone.

For best results use for five minutes, twice a day, for six to eight weeks for maximum benefit.

IMPORTANT SAFEGUARDS

When using ANSR: HALO™ basic safety precautions should always be followed.

READ ALL INSTRUCTIONS CAREFULLY BEFORE USING

DANGER

To reduce the risk of electrical shock:


- Always attach the AC/DC power adapter to the ANSR: HALO™ first, then into an electrical outlet.
- Do not place or store ANSR: HALO™ where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for an ANSR: HALO™ if it has fallen into water while you are recharging the device. Unplug immediately.
- Never operate this device if the AC/DC power adapter or the attached cord has been damaged.

WARNING

To reduce the risk of potentially serious injury:

- Do not aim the red light at your eyes or stare directly at the light source while it is in operation. You may momentarily glance at the light for a second or two. ANSR: HALO™ meets international lamp safety standards for light-emitting diodes.
- Do not allow children to use the ANSR: HALO™ unsupervised.
- Never attempt to open or repair the device or charger. There are no user-serviceable parts inside the ANSR: HALO™. Doing this could result in a serious eye or skin injury or possibly fatal electric shock.
- Use this device only for its intended use as described in this manual. Failure to use and maintain the product in accordance with the instructions outlined in the manual will invalidate the warranty.

SAVE THESE INSTRUCTIONS



MAIN FEATURES & CONTROLS

LENS surface projects red and infrared light and is applied directly to the area to be treated.

POWER BUTTON allows you to turn the device on/off, turn on the red and infrared light, and check the power level of the battery.

BATTERY LEVEL INDICATOR shows three LED lights when fully-charged. Two lights = about 2/3 battery life. One light = 1/3 or less battery left.


POWER JACK PORT is where you plug in the power adapter for charging. The ANSR: HALO will not operate while plugged into an electrical outlet.

GETTING STARTED

- Unpack the ANSR: HALO, the instruction manual, and AC/DC power adapter.
- Read all instructions before using your ANSR: HALO.
- Prior to initial use, charge the device. Use only the AC/DC power adapter provided with the device. Plug the power adapter jack into the bottom of the device and connect to an electrical outlet. Charge ANSR: HALO fully, until all three white indicator lights are illuminated. This will take approximately two hours. Unplug the device once it is fully charged. The ANSR: HALO **will not operate** while plugged into an electrical outlet.

ACTIVATING THE LIGHT

- For best results, wash your face carefully with a quality skin care product before using the light.
 - Activate the light by pressing the power button once.
- The red and infrared LED lights will illuminate along with the battery level indicators. Although not visible (to the




naked eye), the infrared lights will be operating. You will feel ANSR: HALO vibrate while the device is operating.

USING THE LIGHT

- Place the lens in direct contact with your skin and gently move back and forth in a slow, patterned motion over your entire face. Keep the light pressed gently to the skin without pulling or dragging. You may choose to focus on a specific area of concern for the full treatment cycle and then repeat over other areas as necessary. Each treatment cycle is 5 minutes in duration. Complete two treatment cycles a day. The device will automatically shut off after 5 minutes.

CLEANING THE LIGHT

- Use a damp, soft, lint-free cloth and a mild detergent or a disposable, disinfecting wipe to clean the ANSR: HALO lens surface.
- Do not mark, scratch or otherwise damage the lens.





CHARGING THE LIGHT

- Plug the AC/DC power adapter jack into the bottom of the ANSR: HALO and into an electrical outlet. Charge the device fully. Three white lights will illuminate to indicate the device is fully charged. This will take approximately two hours. Unplug the device once it is fully charged.
- To recharge, plug the AC/DC power adapter into the ANSR: HALO and then into the electrical outlet. A single white LED illuminated denotes a charge of up to one-third. Two LEDs indicate about a two-thirds charge. When the device is fully charged, all three indicator lights are lit.

STORING THE LIGHT

- Store your ANSR: HALO in a clean, dry area.



TROUBLESHOOTING

You should reset your ANSR: HALO if one of these things is happening:

- The vibrator is operating without the red light operating
- The device is not going into or coming out of a charge mode
- The battery level indicator lights are behaving strangely

TO RESET THE ANSR: HALO:

Press and hold the power button down until the red light shuts off. This will take approximately 10 seconds. One white indicator light will come on briefly. This will cause the unit to re-enter IDLE mode. Recharge the device.

NOTE: If the reset function does not fix the problem, these symptoms could indicate another hardware issue. In this case, please contact Customer Service at (877) 211-2677.

IF ANSR: HALO WILL NOT TURN ON:

The battery may not be charged. Charge the device for two hours or until all three white indicator lights are illuminated. The product **will not operate** while plugged into an electrical outlet.

IF ANSR: HALO ONLY LASTS A FEW TREATMENTS:

Ensure that you are charging your ANSR: HALO correctly. If the device is fully charged, but still lasts only a few treatments, reset the device and recharge.

IF ANSR: HALO TURNS OFF AFTER FIVE MINUTES:

The five minute treatments are auto-timed and turn off automatically. Push the power button once to restart treatment.

CUSTOMER SERVICE

In North America, call the ANSR: customer service team with any question or to customize your ANSR: skin care regimen. If you are outside North America, contact your local distributor.

WARRANTY

Register your ANSR: HALO at ANSR.com/warranty.

Your ANSR: HALO is warranted to be free from defect in material and workmanship for a period of one year under normal use, after its original purchase date. This warranty extends only to the original retail purchaser with original proof of purchase and only when purchased from an authorized ANSR: retailer.

If the product should become defective within the warranty period, contact ANSR: at (877) 211-2677 for repair or replacement. ANSR: reserves the right to replace a defective product with the most comparable product currently available. Your exclusive remedy for any defective product covered by this limited warranty is limited to the repair or replacement of the defective product. ANSR:'s entire liability for any defective product shall in no event exceed the purchase price for the defective product.

This warranty does not cover products determined by ANSR: to be damaged by any of the following:

- Accident, misuse, abuse or alteration
- Servicing by unauthorized persons
- Use with unauthorized accessories or use other than as instructed
- Connecting to incorrect current and voltage
- Any other conditions beyond our control

ANSR: SHALL NOT BE RESPONSIBLE FOR ANY INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES RESULTING FROM THE USE OF THIS PRODUCT. ALL IMPLIED WARRANTIES INCLUDING BUT NOT LIMITED TO IMPLIED WARRANTIES OF FITNESS AND MERCHANTABILITY ARE LIMITED IN DURATION TO ONE YEAR FROM THE DATE OF ORIGINAL PURCHASE.

This warranty gives you specific legal rights, and you may also have other rights that vary depending on your state, province or other jurisdiction. Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so the above limitations and exclusions may not apply to you.

PATENTS

US Patents Awarded / International Patents Pending

ANSR: HALO

Model Number	ANSR: HALO™
Output Type	Red LED 632nm / Infrared LED 850nm
Battery Type	Rechargeable Nickel Hydride, 1 Cell, AA size
Operating Time	Five minutes per treatment
Operating Temperature	50°F to 104°F (10°C to 40°C)
Operating Humidity	20 to 80% RH non-condensing
Storage Temperature	14°F to 158°F (10°C to 70°C)
Storage Humidity	10 to 95% RH non-condensing
Weight	73.7 Grams
Case Size	78mm x 22mm
Indicators	3 LEDs (Battery level)
Protection	IP44
Equipment Class	Internally Powered

AC/DC Power Adaptor

Model Number	EPS050100-P12-C1
Input Power	110 to 120 VAC, 50 to 60 Hz
Output Power	5 VDC at 1 Amp
Weight	76 Grams
Case Size	63mmL x 26.8mmW x 42mmH
Operating Temperature	32 ⁰ F to 104 ⁰ F (0 ⁰ C to 40 ⁰ C)
Operating Humidity	20 to 80% RH non-condensing
Storage Temperature	14 ⁰ F to 158 ⁰ F (10 ⁰ C to 70 ⁰ C)

Product P/N: _____

ANSR:
PO Box 10087
Portland, OR USA
97296

BE MORE ENLIGHTENED AT [ANSR.COM](https://ansr.com)